MANUAL (p. 2)

ANLEITUNG (S. 8)

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MANUALE (p. 29)

MANUAL DE USO (p. 35)

MANUAL (p. 42)

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KÄYTTÖOHJE (s. 55)

BRUKSANVISNING (s. 61)

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MANUAL DE UTILIZARE (p. 74)

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
## Digital Personal Scale

## Introduction

The digital personal scale is designed to accurately indicate body weight and fat. This scale works on the principle of Bio-electrical Impedance Analysis (BIA). By measuring the electrical resistance and the final calculation of the constant values (age, height, sex, activity level) it is possible to determine the percentage of body fat, body water and muscle percentage in the body.

The scale is able to facilitate 10 persons.

In order to ensure accurate results users should enter the following personal parameters:

- Gender
- Height
- Age
- We recommend that you read this manual before using the scale.
- Keep this manual in a safe place for future reference.

The scale operates with 3x 1.5 V AAA batteries (not included).

## Technical specifications


|                         |                                     |
|-------------------------|-------------------------------------|
| Display:                | Digital LCD display                 |
| Capacity:               | 180 kg/396 lb                       |
| Graduation:             | 0.1 kg/0.2 lb/¼ lb                  |
| Body fat division:      | 0.1%                                |
| Water division:         | 0.1%                                |
| Age range for body fat: | 6–100 years                         |
| Height range:           | 80–220 cm                           |
| Height unit:            | cm/ft-in                            |
| Memory:                 | 10 persons                          |
| Weight units:           | Kilogram (kg)/Stone (st)/Pound (lb) |
| Dimensions:             | 32 x 32 x 2.5 cm                    |
| Weight:                 | 2 kg                                |

## Safety precautions

- To avoid tipping the unit, do not stand on the corners of the platform.

- Do not use the scales if you have a pacemaker or other internal electronic medical device.
- Pregnant women should use the weight function only.
- Do not open the scale in the event of a malfunction. This product should only be opened by a qualified technician.
- Do not expose the scale to high temperatures or humidity.
- Remove the batteries if they will not be used for a long period of time.

## Installation instructions

1. Battery compartment cover
2. Weight units selection button

## Install the batteries

1. Open the battery cover (on the underside of the scale) by releasing the clip and lifting the cover.
2. Insert the batteries.  
☞ Be sure the polarity of the batteries is set properly.
3. Close the cover.

## Select weighing mode


Switch the scale on (ON/OFF button).

Use the Weight units selection button (underside of scale) to select the weight units (kg/lb/st).

## Place the scale

Place the scale on a firm, flat and even surface to ensure accuracy.

## User instructions

1. Display screen
2. Electrodes
3. ON/OFF button
4. SET button
5. UP button
6. DOWN button

## Setting and storing personal data

The scale can store personal data for up to 10 persons.

1. Press the ON/OFF button.  
The personal key number flashes.
2. Press the UP/DOWN buttons to select a personal key (P0 to P9).
3. Press the SET button.  
The gender symbols flash.
4. Press the UP/DOWN buttons to select user gender and condition:
  - Male (♂)
  - Female (♀)
  - Male Athlete (♂ 疾)
  - Female Athlete (♀ 疾)
5. Press the SET button.  
The display shows the height.
6. Use the UP/DOWN buttons to specify user height (user height range may be between 80 cm–220 cm or 2' 7.5"–7' 2.5").

7. Press the SET button.  
The display shows the age (user age range may be between 6 and 100 years).
8. Use the UP/DOWN button to specify user age.
9. Press the SET button.  
The personal data is now set and stored.

## Weighing

1. Press the ON/OFF button.
2. Use the UP/DOWN button to select your personal key number.
3. Step on the scale with bare feet and wait while the scale analyses your weight.
4. On completion, the following screens will be displayed alternately:

1. Weight
2. Fat percentage
3. TBW percentage (Total Body Water)

1. KCAL (daily caloric reference<sup>1</sup>)
2. MUS percentage (muscle)
3. BON (bone weight)

<sup>1</sup> To lose weight, you must consume fewer calories than this.

## Result guidelines – female

### Body fat ratio chart (%)

| Age   | Normal and athlete mode |           |            |           |
|-------|-------------------------|-----------|------------|-----------|
|       | Underweight             | Healthy   | Overweight | Obese     |
| 6–15  | 3.0–16.0                | 16.1–21.0 | 21.1–29.6  | 29.7–50.0 |
| 16–30 | 3.0–16.4                | 16.5–22.0 | 22.1–30.4  | 30.5–50.0 |
| 31–60 | 3.0–17.0                | 17.1–22.4 | 22.5–31.4  | 31.5–50.0 |
| 61–80 | 3.0–17.6                | 17.7–23.0 | 23.1–32.0  | 32.1–50.0 |

Water ratio chart and muscle ratio chart (%)

| Age   | Normal and athlete mode |           |        |                               |
|-------|-------------------------|-----------|--------|-------------------------------|
|       | Low                     | Healthy   | High   | Body Muscle                   |
| 6-15  | < 54.0                  | 54.0-59.0 | > 59.0 | Normal: > 44<br>Athlete: > 46 |
| 16-30 | < 53.0                  | 53.0-58.0 | > 58.0 |                               |
| 31-60 | < 52.0                  | 52.0-57.0 | > 57.0 |                               |
| 61-80 | < 51.0                  | 51.0-56.0 | > 56.0 |                               |

Bone weight chart (kg)

| Weight           | < 45  | 45-60 | > 60  |
|------------------|-------|-------|-------|
| Ref. bone weight | ≈ 3.0 | ≈ 4.2 | ≈ 6.5 |

KCAL Chart

| Age   | Weight (kg) | KCAL (normal mode) | KCAL (athlete mode) |
|-------|-------------|--------------------|---------------------|
| 6-17  | 50          | 1265               | 1445                |
| 18-29 | 55          | 1298               | 1498                |
| 30-49 | 60          | 1302               | 1490                |
| 50-69 | 60          | 1242               | 1332                |
| 70-80 | 50          |                    | 1125                |

Result guidelines – male

Body fat ratio chart (%)

| Age   | Normal and athlete mode |           |            |           |
|-------|-------------------------|-----------|------------|-----------|
|       | Underweight             | Healthy   | Overweight | Obese     |
| 6-15  | 3.0-13.0                | 13.1-17.4 | 17.5-25.6  | 25.7-50.0 |
| 16-30 | 3.0-13.6                | 13.7-18.0 | 18.1-26.4  | 26.5-50.0 |
| 31-60 | 3.0-14.0                | 14.1-18.6 | 18.7-27.0  | 27.1-50.0 |
| 61-80 | 3.0-14.4                | 14.5-19.0 | 19.1-28.0  | 28.1-50.0 |

## Water ratio chart and muscle ratio chart (%)

| Age   | Normal and athlete mode |           |        | Body Muscle                   |
|-------|-------------------------|-----------|--------|-------------------------------|
|       | Low                     | Healthy   | High   |                               |
| 6–15  | < 57.0                  | 57.0–62.0 | > 62.0 | Normal: > 46<br>Athlete: > 48 |
| 16–30 | < 56.0                  | 56.0–61.0 | > 61.0 |                               |
| 31–60 | < 55.0                  | 55.0–60.0 | > 60.0 |                               |
| 61–80 | < 54.0                  | 54.0–59.0 | > 59.0 |                               |

## Bone weight chart (kg)

|        |      |       |      |
|--------|------|-------|------|
| Weight | < 60 | 60–75 | > 75 |
|--------|------|-------|------|

### Maintenance:

Clean the scale with a lightly moistened cloth.  
Do not use cleaning solvents or abrasives.

### Warranty:

Any changes and/or modifications to the product will void the warranty. We cannot accept any liability for damage caused by incorrect use of this product.

### Disclaimer:

Designs and specifications are subject to change without notice. All logos, brands and product names are trademarks or registered trademarks of their respective holders and are hereby recognized as such.

### Disposal:

- This product is designated for separate collection at an appropriate collection point. Do not dispose of this product with household waste.
- For more information, contact the retailer or the local authority responsible for waste management.

This product has been manufactured and supplied in compliance with all relevant regulations and directives, valid for all member states of the European Union. It also complies to all applicable specifications and regulations in the country of sales.

Formal documentation is available upon request. This includes, but is not limited to: Declaration of Conformity (and product identity), Material Safety Data Sheet, product test report.

Please contact our customer service desk for support:

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NEDIS B.V., De Tweeling 28, 5215 MC 's-Hertogenbosch, THE NETHERLANDS (HOLLANDA)